
## "Alles gurgelt!" for all people in Vienna







By showing your barcode, you will receive up to **4 test sets** at any BIPA outlet in Vienna.



Each test set contains **instructions** containing a sample number, a QR code and the link to the **lead-horizon.org/public** WebApp.



By scanning the QR code or entering the URL in a Web browser, you can **register** and activate your sample number. For this, you need to enter your email address and your personal data.



As soon as you are registered, you will receive a **personalised barcode**. With this, you can repeat Step 2 any number of times.



Confirm your identity for the test by holding your ID card or e-card towards the camera.







4





(7)

If you agree to a **video recording** of the test procedure, a **medical lab report** will be issued for you for going to a hair salon or other shops/service providers, if your test is negative. If you skip this step, you will receive a test result, but no officially valid medical report.



Now perform the **test**. You can gargle or rinse your mouth with the saline solution.



Put the **test tube** into the packaging as described and confirm your data. This concludes the test.



Now just hand in the test at any REWE outlet in Vienna (BILLA, BILLA Plus, BIPA and Penny as well as all petrol station shops of BP-Merkur Inside, JET-BILLA Stop & Shop, SHELL-BILLA Unterwegs).



If you hand in the test before 9 a.m., you will receive your result via email within a period of 24 hours. If you have not agreed to the video recording, you can access your test result via the WebApp.



You can access the officially valid medical report via the link you have received by













email.







